

TranSignal™ Human Cytokine Antibody Array 1.0

| | a | b | c | d | e | f | g |
|---|---------------|--------------|-------|---------------|--------------|--------|-----|
| 1 | Eotaxin | GM-CSF | IP-10 | TNF α | Rantes | Leptin | pos |
| 2 | Eotaxin | GM-CSF | IP-10 | TNF α | Rantes | Leptin | pos |
| 3 | IL-1 α | IL-1 β | IL-3 | IL-4 | IL-6 | IL-8 | neg |
| 4 | IL-1 α | IL-1 β | IL-3 | IL-4 | IL-6 | IL-8 | neg |
| 5 | IL-10 | IL-12 | IL-17 | MIP1 α | MIP1 β | MIP-5 | pos |
| 6 | IL-10 | IL-12 | IL-17 | MIP1 α | MIP1 β | MIP-5 | pos |



Isogen Life Science B.V.
 Veldzicht 2A
 3454 PW De Meern
 The Netherlands
 Tel: +31 (0)30 688 0771
 Fax: +31 (0)30 688 8009
 E-mail: info@isogen-lifescience.com

Isogen Life Science N.V./S.A.
 Stokthoekstraat 10A
 9140 Temse
 Belgium
 Tel: +32 (0)9 351 48 79
 Fax: +32 (0)9 277 93 03
 E-mail: info@isogen-lifescience.com

Isogen Life Science S.L.
 Poligono Les Guixérés-Edif. BCIN
 C/ Marcus Porcius 1
 08915 Badalona Spain
 Tel: +34 934 64 80 07
 Fax: +34 934 64 80 20
 E-mail: clientes@isogen-lifescience.com